

## EYEDAPTIC TECH NOTE: SETTING FAVORITES

The Eyedaptic® app allows the user to set up to 3 favorites and switch between them at any time.


These favorite settings can be used to save “modes” for repeated tasks. For example, if the user watches TV in the same spot, with the same lighting conditions, a “TV watching mode” could be set up as one of the favorites by saving the desired magnification and brightness settings as Favorite #1. Then, every time Favorite #1 is selected, these settings will automatically be adjusted.

### How to set Favorites:


1. Adjust the features to the desired settings
  - a. For example:
    - i. Set magnification to 5x
    - ii. Set brightness to 85%
2. On the Eyedaptic® app, press “Save as” underneath the Favorites options




3. Then press the desired favorite number, ie. Favorite #1



4. Then press “Save”



5. The features are now saved under this favorite number, so every time that number is selected, these features will be automatically adjusted
- a. When a favorite is activated, the number button will be black
  - b. When a favorite is not activated, the number button will be white



\*Note: favorites can always be changed and overridden by completing steps 1-5 again

\*Note: when a feature is adjusted, ie. magnification is increased, the favorite will “turn off” – it can be turned back on again by simply pressing the favorite number